

Kroehling Advanced Materials Foundry

Bruker Q4 Tasman Advanced CCD-Based Optical Emission Spectrometer

Spectrometer Elements

Iron & Steel Alloys: C, Si, Mn, P, S, Cr, Mo, Ni, Cu, Al, As, B, Bi, Ca, Ce, Co, Nb, Pb, Sb, Sn, Ta, La, Ti, V, W, Zr

Aluminum Alloys: Si, Fe, Cu, Mn, Mg, Cr, Ni, Zn, Ti, Ag, B, Ba, Be, Bi, Ca, Cd, Co, Ga, In, Li, Mo, Ni, P, Pb, Sn, Sr, V, Zr, Sb, Hg

Copper Alloys: Zn, Pb, Sn, P, Mn, Fe, Ni, Si, Mg, Cr, Al, S, As, Be, Ag, Co, Bi, Cd, Sb, Zr, Ti, Au, C, Nb, Se, Te